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VIVIT IF Adipose-Rich Tissue Clearing Kit (Ionic Liquid-Based)

Product information

The VIVIT Tissue Clearing Kit series, leveraging innovative ionic liquid-based clearing technology, establishes a third mechanism for tissue refractive index matching, achieving both "high transparency" and "minimal deformation". It uniquely offers enhanced fluorescence and eliminates ice crystal formation at low temperatures (Cell, 2025), providing an innovative solution for 3D imaging of thick ex vivo tissues. The kit contains core components including delipidation reagent, overionization reagents, and ionic liquids matching reagent. These components transform tissues into an ionic glass-like state, achieving optical transparency, enabling highly transparent tissue clearing. Critically, tissues processed this way can be stored long-term at low temperatures (-80℃) without fluorescent signal decay. VIVIT-processed tissues are compatible with various 3D fluorescence imaging methods. Featuring user-friendly operation and stable outcomes, it provides a reliable tool for tissue analysis in scientific research. For immunofluorescence-stained adipose-rich tissues, such as viscera, the VIVIT IF Adipose-rich Tissue Clearing Kit is recommended.

Reagents and Storage Condition

Reagent	Size	Notes	Storage
Antibody cross-linking reagent	1 mL*3	For cross-linking and fixation of antibodies to tissues after	-80℃ or -20℃
	1 mL*1		
Delipidation reagent A2 (DIDC-C)	30 mL	For delipidation of tissue	Store at room temperature.
	10 mL		
Delipidation reagent B	30 mL	For gradient delipidation of adipose-rich tissue	Store at room temperature.
	10 mL		
Overionization reagents B (SDS)	30 mL	For overionization of adipose-rich tissue	Store at room temperature.
	10 mL		
Ionic liquids matching reagent (RIL-C)	30 mL	Ionic liquid solution for rinsing purposes	Optimal temperature: -80℃ Do not store at 4℃.
	10 mL		
Ionic liquids matching reagent (IL5-C)	30 mL	Ionic liquids matching reagent	
	10 mL		
Ionic liquids matching reagent (IL2.5-C)	30 mL	Ionic liquids matching reagent	
	10 mL		

Required material not supplied

0.01 M PBS

Tissue Processing Guidelines

Prior to tissue clearing, tissue fixation should be performed either by passive immersion or by intracardiac perfusion. Perfusion fixation with 4% PFA is recommended. After fixation, thoroughly rinse the tissue with 0.01 M PBS to remove residual PFA. The kit is designed for ex vivo tissue clearing following immunofluorescence staining. The immunofluorescence staining can be performed using either conventional immunofluorescence staining protocols or thick tissue immunofluorescence staining methods (such as whole-mount, SWITCH, eFLASH, etc.).The tissue clearing process involves steps such as delipidation, overionization, and ionic liquids matching. If issues arise during processing, refer to the [Precautions] section.

Specific steps (suitable for thick tissues of 1mm thickness; reduce processing time for thinner tissues and increase for thicker tissues as appropriate):

- Antibody Cross-linking: Transfer the tissue that has completed immunofluorescence staining into the antibody cross-linking reagent (diluted 10-fold with 0.01 M PBS solution). Incubate at 4℃ for 8-12 hours.
- Tissue Delipidation: Transfer the tissue to delipidation solution A2 and incubate on a shaker at 37℃ for 2-4 hours. Then transfer the tissue to a mixture of 20% delipidation solution A1 and 80% delipidation solution B, and incubate at room temperature for 2-4 hours. Subsequently, transfer the tissue to delipidation solution B and incubate at room

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
temperature for 2-4 hours. Finally, transfer the tissue to delipidation solution A1 and incubate at 37℃ for 2-4 hours. The processing time may be appropriately extended based on the observed effect.

- Tissue Overionization: Discard the delipidation reagent A2 in the container. Add overionization reagent B. Incubate on a shaker at room temperature for 4 hours.
- Rinsing: Discard the reagent in the container and replace it with ionic liquid rinsing reagent RIL-C pre-equilibrated to room temperature. Invert to mix immediately after adding. Perform an ice-bath rinse on a shaker for 1 hour; additionally, inverting and manually mixing several times during this step is recommended.
- Refractive Index Matching: Discard the reagent from the container and replace it with ionic liquid matching reagent IL5-C. Invert to mix immediately after adding. Incubate on a shaker at room temperature for 2-4 hours; additionally, inverting and manually mixing several times during incubation is recommended. Subsequently, discard the reagent and replace it with ionic liquid matching reagent IL2.5-C. Immediately after addition, invert the container to mix thoroughly. Incubate on a shaker at room temperature for 4-6 hours, with additional manual inversions recommended during processing, until the tissue achieves complete transparency.

Precautions

- This kit is intended solely for the optical clearing of ex vivo tissues and must not be used on in vivo tissues.
- This kit has a one-year shelf life; the manufacturing date can be found on the product label.
- All reagents stored at low temperatures must be equilibrated to room temperature before unsealing and use. Reagents may be stored at room temperature if used within one week.
- Upon arrival, it is recommended to aliquot and store the antibody crosslinker. This reagent is prone to degradation upon exposure to moisture. After dilution with PBS, use immediately. The diluted solution will be severely compromised if stored at room temperature for more than 2 hours.
- IL2.5-C and IL5-C may crystallize when stored at 4℃ or -20℃. Storage at -80℃ prevents crystallization of these reagents.
- If crystallization occurs in IL5-C or IL2.5-C, tightly seal the reagent vessel and incubate in a 50℃ water bath or metal bath until crystals fully dissolve (approximately 30 minutes). The reagent can be used normally after complete dissolution.
- Processing duration for each step may be extended based on experimental outcomes. Generally, increase all step durations for thicker tissues or samples subjected to extended fixation.
- To ensure optimal tissue clearing quality, use an orbital shaker or intermittent manual mixing after tissue immersion in reagents.
- Recommended reagent volume: ~3 mL per tissue sample (1 cm² × 1 mm thickness). Insufficient volumes may compromise clearing efficacy.
- VIVIT-processed tissues may be stored long-term at -80℃.
- Pausing points may be set at any step except during RIL-C rinsing. Extending other steps (excluding RIL-C rinsing) has minimal impact on results. Overnight processing is recommended for delipidation and IL2.5-C refractive index matching steps.
- DIDC-C effectively removes oil-based pen marks. Always close lids tightly during use to prevent leakage.
- This product is restricted to scientific research by qualified professionals only. Not for clinical diagnosis or treatment.
- For your safety and health, please wear disposable gloves and a mask when using this reagent.
- Keep all reagents away from fire and heat sources, and store them in a cool, dry place below 25℃, avoiding direct sunlight.
- All waste liquids from tissue processing should be collected separately in dedicated containers. Do not mix them with strong oxidizing agents, and avoid pouring waste liquids directly into drains.
- First-time users: Follow the specified volumes/times in the manual and purchase optional reagents as needed. Users with special requirements should contact technical support via our official website.

效果示意 Effect Illustration

参考文献 Reference

VIVIT: Resolving trans-scale volumetric biological architectures via ionic glassy tissue
<https://doi.org/10.1016/j.cell.2025.07.023>